

“Neck Extender” for **Priefert**® Model 91 Headgate Mounting and Installation Instructions (HGNE)



Priefert's Neck Extender for their Model 91 Headgate is designed to provide improved access to the neck area when vaccinating and to help restrict movement of the cow's head aiding procedures such as tagging, dehorning, eye medicating, etc.



CAUTION! FOR SAFEST INSTALLATION...

Headgate yokes should be blocked or chained in the OPEN position to prevent unexpected closure during installation. This measure also allows for easy access to key parts of the Model 91 Headgate during assembly of neck extender.

Mounting Extender to Existing Headgate:

Tools Required: 9/16" Wrench



Step 1: Locate recommended position for right and left Extenders; safety shields at the bottom as seen above. Position top of extender by installing U-Bolt around top of Yoke, just below roller box (Fig. 1). Secure Extender with two (2) U-Bolt nuts. Finger tighten only.

Step 2: Slide Extender upward and slip U-bolt around bottom of Yoke just above slide collar (Fig. 2). Attach Extender with two (2) U-Bolt nuts. Finger tighten only.

Step 3: To insure proper alignment, evenly tighten each nut on the top of the Extender until snug. (Fig. 3) **DO NOT OVERTIGHTEN.**

Step 4: Evenly tighten bottom nuts in the same manner. (Fig.4)

When proper alignment is achieved, tighten all nuts again.

Step 5: Repeat Steps 1-4 for other side.

Benefits of Priefert's Accessory Neck Extender for the Model 91 Headgate:

- Provides 5" of easy neck access.
- Helps hold the in an extended position to make vaccinations easier and other procedures around the head safer.
- Attaches easily with included U-Bolts.
- Architectural Grade Powder Coat Finish.
- Designed to mount exclusively to Priefert's Model 91 Headgate.
- Can be purchased as an add-on to a new Model 91 headgate or chute. Easily added to an existing chute if chute has Priefert Model 91 Headgate.



Fig. 1
Top U-Bolt



Fig. 2
Bottom U-Bolt



Fig. 3
Tightening
Top Nuts



Fig. 4
Tightening
Bottom Nuts

2630 S. Jefferson
P.O. Box 1540
Mt. Pleasant, TX 75456-1540

Priefert®
Rodeo & Ranch Equipment

1-800-527-8616
903-572-1741 903-572-2798 Fax
sales@priefert.com